

To protect the privacy of
customers or suppliers of
Salamander Jewelry.
We erased their address from
this document.

Document No. : TR-EO-002181116

Date : 18/11/2016

Page : 2 / 4

PRODUCT : ORDER NO. 2166123

SAMPLE SLIP PARAMETER

EXHAUST:04H 00MIN TOTAL :11H 56MIN
PRETT :00H 50MIN STERI :07H 00MIN
ACTUAL TIME

=====

05:28:20 STERILIZATION SUCCESS

END 05:28:12 058,40 -004kPa
05:28:12 058,40 -004kPa
01:01:24 058,60 -015kPa
EXHAUST 01:22:12 058,50 000kPa
01:21:42 058,20 000kPa
01:11:42 058,30 000kPa
01:01:42 058,30 000kPa
00:51:42 058,30 002kPa
00:41:42 058,20 000kPa
00:31:42 058,30 001kPa
00:21:42 058,30 000kPa
00:11:42 058,20 001kPa
00:01:42 058,40 000kPa
23:51:42 058,90 001kPa
23:41:42 058,40 001kPa
23:31:42 058,30 000kPa
23:21:42 058,40 000kPa
23:11:42 058,40 000kPa
23:01:42 058,40 000kPa
22:51:42 059,00 001kPa
22:41:42 058,30 000kPa
22:31:42 058,50 000kPa
22:21:42 058,70 003kPa
22:11:42 058,90 000kPa
22:01:42 058,90 000kPa
21:51:42 058,90 000kPa
21:41:42 059,00 000kPa

21:31:42 059,00 000kPa
21:21:42 058,90 001kPa
21:11:42 059,00 000kPa
21:01:42 059,00 001kPa
20:51:42 059,00 001kPa
20:41:42 058,90 001kPa
20:31:42 059,00 057kPa
20:21:42 059,00 060kPa
20:11:42 059,10 061kPa
20:01:42 059,10 062kPa
19:51:42 059,20 067kPa
19:41:42 059,30 075kPa
19:31:42 059,20 087kPa
19:21:42 059,10 100kPa
19:11:42 059,20 117kPa
19:01:42 059,10 146kPa
18:51:42 059,10 219kPa
18:41:42 059,30 311kPa
18:31:42 059,30 311kPa
STERI 18:22:33 059,50 224kPa
18:21:42 059,60 010kPa
ADD EO 18:17:39 059,70 -087kPa
18:11:42 059,10 -017kPa
HUMID 18:07:57 059,20 -067kPa
PRETT 18:01:42 059,20 001kPa
PRICAT 17:32:46 055,90 000kPa
STATE TIME TEMP PRE

=====

RUN COUNT : 0965
EO PRE :060kPa VENT TIME:04H 00M
STERI TEMP: 058,00 STERI TIME:07H 00M
BRO TIP 17:32:46 055,90 000kPa
STATE TIME: 2016-11-17 17:32:38